



Pure Comfort

3-Ball Breath Exerciser

- Wide flow rate from 600-1200cc/Sec.
- Compact and portable
- Break-resistant plastic
- Color coded flow rate visualization
- Improves airflow and optimizes lung capacity.
- Strengthens primary respiratory muscles.
- Measures the volume of air inspired and expired by the lungs.
- Allows visual calibration and estimation of the user's progress.
- Single Patient Use
- Latex Free
- Non-Sterile



ID# 60006-0377-83



Distributed by:
Home Aide Diagnostics, Inc.
1072 S. Powerline Rd.
Deerfield Beach, FL 33442
1-800-915-0116
www.homeaide.us