

Pro Comfort

PULSE OXIMETER



NON-INVASIVE MEASUREMENT OF ARTERIAL BLOOD OXYGEN SATURATION AND PULSE RATE



PRBPM
PARAMETER


AUTO
SHUT-OFF
LOW POWER
INDICATOR

%SPO2
HEMOGLOBIN
SATURATION

NDC # 50632-0007-53

- Lightweight for carrying and Easy-To-Use.
- Manually adjust the direction of interface.
- LED display, simultaneous display for testing value and bargraph.
- Low Battery voltage indicator.
- Automatically standby or sleep.

OLED display modes



1. Display Type: OLED display

2. SpO2:

Measurement range: 70%-99%

Accuracy: $\pm 2\%$ on the stage of 70%-

99%, Unspecified (70%)

Resolution: $\pm 1\%$

3. PR:

Measurement range: 30BPM-240BPM

Accuracy: $\pm 1\text{BPM}$ or $\pm 1\%$

4. Power Supply: 2 AAA 1.5V batteries

5. Dimension: 60*35*35mm



Home Aide Diagnostics
1072 S. Powerline Rd.

Deerfield Beach, FL 33442

www.homeaide.us | 1-800-915-0116