




Pro Health Track

Bluetooth Blood Pressure Monitor



Model: DBP-1332b



2 Users



Auto power-off
Low battery detection

Upper Arm Cuff 8.66-14.17in (22-36cm)



Oscillometry results:

- Systolic (60~280) mmHg
- Diastolic (30~200) mmHg
- Pulse (30~180) BPM
- mmHg or kPa
- Hypertension
- Hypotension
- 120 Memories

Name	ID	UPC
Pro Health Track: Bluetooth Blood Pressure Monitor	50027-0494-45	850027494456



Home Aide Diagnostics, Inc.
 1072 S. Powerline Rd.
 Deerfield Beach, FL 33442
 1-800-915-0116
 www.homeaide.us