

# TrueComfort

## SAFETY PEN NEEDLES

- ☀ 31 OR 32 GAUGE
- ☀ 100 PER BOX
- ☀ SINGLE USE
- ☀ SAFETY CAP
- ☀ STERILE
- ☀ UNIVERSAL SIZE



Clear safety shield automatically locks after 1 use to protect the patient

Ergonomic base with ribbing for comfortable grip



Visual indication when the needle is permanently locked

Compatible with most pen needle devices

Product Name	UPC	ID	Box QTY	Case QTY	AWP
True Comfort Safety Pen Needles 4MM 32G	850027494449	50027-0494-44	100/box	20	\$144.87
True Comfort Safety Pen Needles 5MM 31G	850027494432	50027-0494-43	100/box	20	\$144.87
True Comfort Safety Pen Needles 6MM 31G	850027494425	50027-0494-42	100/box	20	\$144.87



Home Aide Diagnostics, Inc.  
 1072 S. Powerline Rd.  
 Deerfield Beach, FL 33442  
 954-794-0206  
[www.homeaide.us](http://www.homeaide.us)